BEAUTIFLY

Glossy

Blush

series

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Rated voltage:	5V/1A
Rated power:	4W
Battery Capacity:	600mAh
Charging time:	3 hours
Massage head diameter:	35mm

Description of the device [page 2, f1]

1. Ultrasonic head
2. Photon lights
3. Battery charge level
4. Phototherapy: Red/Blue
5. Operating modes
6. On/off button
7. Mode change button

4 modes – ways of using:

Note: Start your care with cleansing (Cleanse mode), use a facial cleansing cosmetic. Thorough cleansing of pores will involve deep removal of impurities, and the dirt will come to the surface of the skin. Wash the skin surface with make-up remover. Set each mode for 4 minutes of use time. Within 1 minute of inactivity, the device will automatically turn off.

The device works on the induction principle, which means that it starts working after contact with the skin.

Using Cleanse Mode (using a cotton pad)

In this mode you can choose from 3 levels of intensity, heat and vibration. Thanks to this, you can effectively cleanse pores of impurities, brighten the skin and improve its color.

1. Remove the hook on the cotton pad
 2. Attach a cotton pad
 3. Pour makeup remover onto the pad
 4. Press and hold the button for 3 seconds to turn on the device
 5. Place your fingers around the sensor strips
 6. Cleanse your face in circular motions
- CARE TIME: 4 minutes

Applying Rejuvenate Mode (no need to use a cotton pad)

The mode uses red light and heat and vibration massage, as well as blue light, heat and vibration massage. It has a strong firming effect, stimulating the skin to intensively produce collagen and elastin. Effectively eliminates wrinkles and furrows, moisturizes, leaving the effect of rejuvenated skin.

1. Apply a serum or mask to your skin
 2. Press the power button for 3 seconds to turn on the device and select Rejuvenate mode
 3. While holding the device, touch the sensor strips on both sides of the device. Perform the treatment by pulling the device from the bottom of the face upwards
- CARE TIME: 4 minutes

Using Lifting Mode (no need to use a cotton pad)

This mode effectively stimulates and relaxes the skin and helps deeply moisturize it. Red light combined with vibration and heat greatly supports facial remodeling, allowing it to regain a beautiful, youthful oval.

1. Press and hold the power button for 3 seconds to turn on the device. Select Lift mode
 2. While holding the device, touch the sensor strips on both sides of the device. Drag the device from the bottom of your face to the top
- CARE TIME: 4 minutes

Using Eye care mode

The mode uses vibration, heat and optional red light. After facial skin treatments, we move on to the eye area and perform precise care of the delicate skin around the eyes. The skin around the eyes requires special care, and appropriate treatments will help you eliminate facial wrinkles, swelling and dark circles under the eyes.

CARE TIME: 4 minutes

Instructions for use

Using the power adapter

1. When the battery symbol appears on the controller display this means the battery is low. Insert the plug of the power adapter into the DC socket of the controller, then insert the plug of the controller into the power socket, the battery indicator on the LCD will flash.
2. The battery life of the rechargeable battery is equal to or greater than 500 charge-discharge cycles under normal conditions, which allows it to be used for more than three years (once a day).

Precautions

1. Do not leave the device connected to the power supply for a long time.
2. Do not touch it with wet hands, it may result in electric shock or short circuit.
3. Do not place the device in places with high humidity, as it may result in electric shock, short circuit or insulation damage.
4. Do not leave the device within the reach of children.
5. Children are prohibited from using the device.
6. It is prohibited to use it in places where flammable materials are stored.
7. Do not clean the device with corrosive substances.
8. Use only the included power supply.
9. Before connecting the power supply, make sure the device is dry.
10. Do not use the device on sunburned skin.
11. Discontinue use in case of unwanted side effects.

Contraindications

- pregnancy
- metal implants, pacemaker
- metal objects (e.g. jewelry)
- tumor
- fever
- disease states, weakness
- disruption of the continuity of the epidermis
- previously used fillers, botox (6 weeks must have passed since the treatment)
- viral, bacterial and fungal skin diseases
- acute forms of rosacea and phlegmon
- active herpes
- epilepsy
- migraines
- untreated and unregulated thyroid diseases.

More information about the product is available on the manufacturer's website at: [page 2, f2]

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[EN] Rules for handling electrical and electronic equipment: Electrical and electronic equipment as well as batteries and accumulators must not be combined with other household waste. This is indicated by the crossed-out wheeled bin symbol on the packaging. Dispose of waste electrical and electronic equipment as well as batteries and accumulators at your local separate collection point or to your dealer. Electrical and electronic equipment and batteries and accumulators may contain hazardous substances, mixtures and components that are harmful to the environment and human health. Incorrect handling or damage to them may lead to damage to health or environmental pollution when the equipment is further disposed of or recycled. The device was placed on the market after August 13, 2005

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[EN] Summary declaration of conformity: The manufacturer hereby declares under its sole responsibility that the above product complies with Directive 2014/30/EU of the European Parliament and of the Council of 26 February 2014 on the harmonisation of the laws of the Member States relating to electromagnetic compatibility and Directive 2011/65/EU of the European Parliament and of the Council of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment, together with Commission Delegated Directive (EU) 2015/863 of 31 March 2015 amending Annex II to Directive 2011/65/EU of the European Parliament and of the Council as regards the list of restricted substances.

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Importer: Beautifly Sp. z o.o.

Plac Bankowy 2, 00-095 Warszawa, Polska

Model: Glossy Blush